

流式细胞仪	激光器	检测器
BD LSRFortessa承担所有分析实验	488 nm	695/40
		530/30
	405 nm	780/60
		710/50
		670/30
		610/20
		525/50
		450/50
	561 nm	670/30
		610/20
582/15		
BD FACSAria™ III 用于动植物细胞分选	488 nm	695/40
		530/30
	405 nm	510/50
		450/40
	561 nm	780/60
		670/14
		610/20
	633 nm	582/15
780/60		
660/20		
BD FACSAria™ SORP用于微生物细胞分选	488 nm	695/40
		530/30
	355 nm	675/50
		450/50
	405 nm	525/50
		450/50
	561 nm	780/60
		710/50
		610/20
	633 nm	582/15
780/60		
670/30		
Beckman MoFlo XDP用于哺乳动物细胞分选	488 nm	529/28
		575/25
		620/29
		670/30
		785/62
	355 nm	457/50
		620/29
	642 nm	670/30
		670/30
		720/13